



Rumen Magnet with Plastic Cage



Weiyou® Rumen Magnet with plastic cage manufacturer and supplier in China. We have been specialized in veterinary products for many years. Weiyou® Rumen Magnet with plastic cage have a low-price advantage and cover most of the European and American markets. Weiyou® look forward to becoming your long-term partner in China.

Product Description

Weiyou® Rumen Magnet with plastic cage is high quality Ferrite and ABS plastic. Rumen Magnet with plastic cage is mainly used for cattle and cow, and is rod magnet with rounded ends, which are popular with dairy farmers and veterinarians to help prevent hardware disease. The magnets have the correct size, shape and weight to offer maximum protection and cure against traumatic indigestion in cattle. Color of plastic is available, and OEM packaging is acceptable.

Specification

size	100x25mm,100x37mm,
Material	Ferrite and ABS
temperature	-30~80
Usage	cattle
Shelf life	3 years

Features

Used worldwide as an effective remedy for hardware disease.

Classic northern hemisphere style with spaces for accumulating hardware within the plastic cage.

Color of plastic is available

Packaging: 10pcs with one neuter box, 200pcs per carton

OEM packaging is available

Size: 100X37CM, 100X25CM







www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone:+86-13567833512
E-mail:_dario@nbweiyou.com



www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone:+86-13567833512
E-mail:_dario@nbweiyou.com



